

Andreasen pipette

BS 812

GENERAL DESCRIPTION

The Andreasen pipette is used to extract precise quantities of suspension ready for analysis. The Pipette is positioned on the pipette stand that let the vertical movement to be conveniently adjusted without disturbing the suspension under test.

ORDERING INFO

47-D0439/B

Andreasen pipette, 25 ml capacity

ACCESSORIES

22-T0062/2A

Pipette stand with millimetre scale.


Weight: 10kg.



Andreasen pipette (47-D0439/B)
and stand (22-T0062/2A)



Andreasen pipette stand 22-
T0062/2A



Andreasen pipette millimetre scale