

DATA SHEET

C NTROLS

Kumagawa extractor

EN 12697-1

GENERAL DESCRIPTION

Used for the quantitative determination of bitumen in hot-mixed paving mixtures and pavement samples. It consists of a balloon, cooling unit, Dean-Stark stopcock, electric heating mantle with regulator and fittings. Two models available: 1 liter and 2 liter capacity. The filtering cartridges are not part of the apparatus and have to be ordered separately. See accessories.

The solvent used can be profitably recovered with the Solvent recovery unit.

• Weight approx.: 15 kg

ORDERING INFO

75-B0018

Kumagawa extraction apparatus, 1 liter capacity. 230 V, 50-60 Hz, 1 ph

75-B0018/Z

Kumagawa extraction apparatus, 1 liter capacity. 110 V, 60 Hz, 1 ph

75-B0018/A

Kumagawa extraction apparatus, 2 liters capacity. 230 V, 50-60 Hz, 1 ph

75-B0018/AZ

Kumagawa extraction apparatus, 2 liters capacity. 110 V, 60 Hz, 1 ph

ACCESSORIES

75-B0018/1

Kumagawa cartridges 58 mm dia. x 170 mm, for 75-B0018 - 1 liter capacity extractor. Pack of 25

75-B0018/A1

Kumagawa cartridges 80 mm dia. x 200 mm, for 75-B0018/A - 2 liters capacity extractor. Pack of 25

75-B0018/A2

Kumagawa cartridges 80 mm dia. x 250 mm, for 75-B0018/A - 2 liters capacity extractor. . Pack of 25



75-B0018